

201-C Electrode



Main Features

- Ideal for measurements in all common aqueous solutions (scientific research, environmental engineering, quality control)
- High-quality pH glass electrode
- Removable storage cap retains moisture
- KCl gel electrolyte; no refilling required
- Not suitable for solutions with low ionic concentration (e.g. deionized/distilled water), strongly alkaline solutions (pH>12), highly corrosive liquids, or long-term measurements at high temperature (>60°C)
- Compatible with all pH meters with a BNC connector

Technical Parameters

Measuring range	(0 ~ 14) pH	Connector	BNC
Shaft material	PC (polycarbonate)	Cable	Ø3×1 m
Junction	Single ceramic	Electrode dimensions	(Ø12×160) mm
Reference system	Ag/AgCl	Temperature sensor	No
Electrolyte	Gel KCl	Recommended operating temperature	(0 ~ 80)°C