

2401-C Electrode



Main Features

- Specially developed for high-precision laboratory measurements in the 0 to 200 mS/cm range
- Cavity structure and platinum chip sensor design improve measurement accuracy and stability
- Glass body resists chemical corrosion
- No integrated temperature sensor
- ATC requires a separate temperature probe or manual temperature entry
- Compatible with all conductivity meters with a BNC connector

Technical Parameters

Measuring range	20 to 2000 mS/cm,	Cell constant	K = 1,0
Shaft material	lead-free glass	Connector	BNC
Temperature sensor	No	Cable	1 m
Recommended operating temperature	(0 ~ 80)°C	Electrode dimensions	(Ø12×145) mm