

DJS-0.1-F Electrode



Main Features

- Integrated temperature sensor, especially suited to measuring the conductivity of purified water and ultrapure water
- Flow-through module eliminates the influence of atmospheric CO₂
- Highly sensitive platinum chip sensor design for high accuracy at low conductivity (0 to 200 $\mu\text{S}/\text{cm}$)
- Integrated temperature sensor enables automatic temperature compensation
- Compatible with all conductivity meters with a BNC connector

Technical Parameters

Measuring range	0 to 200 $\mu\text{S}/\text{cm}$	Cell constant	K = 0,1
Shaft material	lead-free glass	Connector	BNC/RCA
Temperature sensor	30 k Ω thermistor	Cable	1 m
Electrode	$\varnothing 7 \times 18$ mm platinum	Electrode dimensions	($\varnothing 12 \times 155$) mm