

EC60-DE Replacement Electrode



Main Features

- Compatible with the EC60-Z, PC60-Z, EC60, and PC60 meters
- Simultaneous measurement of conductivity, TDS, salinity, resistivity, and temperature
- Electrode with a durable platinum black coating – precise and robust

Technical Parameters

Conductivity	Measuring range	0 to 20 mS/cm
	Cell constant	K=1
	Sensor type	Platinum black sensor
Temperature	Temperature sensor	30 kΩ thermistor
	Temperature compensation	0 ~ 50°C (32 ~ 122°F)
Other	Body	Polycarbonate
	Dimensions	75×35×35 mm
	Weight	38 g