

LabSen 551 Electrode



Main Features

- Patented Swiss sensor technology and materials
- Pointed electrode for direct measurement of soil, growing media, fruit, vegetables, etc.
- Pointed glass membrane delivers stable readings in 15–30 seconds, up to 10 times more robust than conventional electrodes
- Polymer electrolyte with pinhole junction structure prevents clogging; no refilling required
- PVC body is sturdy and resistant to fouling
- Long-life reference system (glass tube, AgCl, silver reference wire) improves stability and service life
- Compatible with all pH meters with a BNC connector

Technical Parameters

Measuring range	(0 ~ 14) pH	Electrolyte	Polymer
Recommended operating temperature	(0 ~ 80)°C	Storage solution	3M KCl
Shaft material	PVC plastic	Membrane impedance	<250 MΩ
Glass membrane type	S	Electrode dimensions	Electrode (Ø15×90) mm Sensor tip (Ø6×25) mm
Reference system	Long-life	Connector	BNC
Junction	Ceramic + single pore	Cable	Ø3 mm×1 m