

Magnetic Stirrer with Electrode Holder 601



Main Features

- Compact and rugged magnetic stirrer
- Continuously variable speed control
- Electrode holder for mounting measuring electrodes
- includes two stir bars (25 mm and 30 mm)
- Speed range: 0 – 2300 rpm
- Work plate diameter: 120 mm
- Electrode clamp dimensions: $\varnothing 16$ mm $\times$ 13 mm
- Electrode clamp height: 0 to 99 mm
- Max. stirring volume: 3000 ml
- Power adapter: DC 9V
- Unit dimensions: $\varnothing$ 120 mm $\times$ 49 mm
- Weight: 500 g