

Heated Magnetic Stirrer 802



Main Features

- Easily adjustable speed (0 to 2300 rpm, max. stirring volume: 1000ml)
- Heats the stirring surface to 120°C (248°F) within 5 minutes
- Overheat protection stops the heating when the stirring surface exceeds 150°C (302°F).
- The heat-resistant magnetic steel effectively prevents stir bar hopping.
- Well suited for stirring low-viscosity solutions
- Work plate diameter: 120 mm
- Hot plate temperature: 120°C to 150°C
- Dimensions: Ø 120 mm × 53 mm
- Weight: 500 g