

ORP60-DE Replacement Electrode



Main Features

- Compatible with the ORP60 and ORP60-Z meters
- Platinum ORP sensor
- Double junction structure
- Integrated temperature sensor
- Suitable for ORP testing of common aqueous solutions

Technical Parameters

ORP	Measuring range	-1000 ~ 1000 mV
	Sensor type	Platinum
	Junction	Double ceramic
	Electrolyte	Gel KCl
	Reference	Ag/AgCl
Temperature	Temperature sensor	30 k Ω thermistor
	Operating temp.	0 ~ 50°C (32 ~ 122°F)
Other	Body	Polycarbonate
	Dimensions	75×35×35 mm
	Weight	40 g