

PC60-E Replacement Electrode



Main Features

- Compatible with the PC60 meter
- Simultaneous measurement of pH, conductivity, and temperature
- Measurement of pH, conductivity (EC), TDS, salinity, and temperature
- Unique BPB (Brush-resistant Platinum Black) sensor allows a wide conductivity measuring range with high accuracy

Technical Parameters

pH	Measuring range	0 to 14
	Reference	Ag/AgCl
	Junction	Ceramic
	Sensor type	Lithium glass
	Electrolyte	3M KCl
Conductivity	Measuring range	0 to 20 mS/cm
	Cell constant	K=1
	Sensor type	Platinum black sensor
Temperature	Measuring range	0 to 50 °C
	Temperature sensor	30 kΩ thermistor
Other	Body	Polycarbonate
	Dimensions	75×35×35 mm
	Weight	40 g