

PC60 Premium 5-in-1 Multi-Parameter Pocket Meter



Main Features

- **5-in-1 Meter** for pH, Conductivity, TDS, Salinity, and Temperature
- **Simultaneously tests** pH / EC / TDS / Salinity / Temperature quickly, accurately, and reliably
- Replaceable multi-parameter electrode
- **Large LCD display** with 2 backlight colors (blue for measurement mode and green for calibration mode)
- Easy calibration with **automatic buffer recognition**
- **Quick switching** between pH, conductivity, TDS, and salinity measurement modes
- **Automatic Temperature compensation** (ATC) (range 0-50°C)
- **Complete kit** incl. calibration solutions, storage solution, batteries, and carrying strap in a convenient carrying case



Technical Parameters

pH	Measuring range	-2,00 to 16,00pH
	Resolution	0,01pH
	Accuracy	± 0,01pH (±1 digit)
	Calibration Points	1 to 3 points
	Temp.-Compensation	0 to 50 °C (automatic)
Conductivity	Measuring range	0 – 200,0 µS, 200 – 2000 µS, 2 – 20,00 mS/cm
	Resolution	0,1/1 µS, 0,01 mS/cm
	Accuracy	±1 % F.S
	Calibration Points	1 to 3 points
TDS	Measuring range	0,1 ppm to 10,00 ppt
	TDS-Factor	0,40 to 1,00
Salinity	Measuring range	0 to 10,00 ppt
Temperature	Measuring range	0 to 50 °C
	Accuracy	±0,5 °C
Other	Power supply	4 × AAA Batteries
	Operating time	Up to 1000 hours
	IP rating	IP67
	Dimensions	40×40×178 mm (with case 255×210×50 mm)
	Weight	133 g (with case 800 g)
	Auto. Power-Off	In 8 minutes without operation
	Various display modes	Yes
	Display	2-Color Backlit LCD

What's in the Kit

PC60 Meter	1
pH4.00 Calibration solution (50 ml)	1
pH7.00 Calibration solution (50ml)	1
1413µS/cm Calibration solution (50ml)	1
12.88mS/cm Calibration solution (50ml)	1
Calibration vial	4
3M KCl Electrode Storage Solution	1
AAA-Batteries	4
Manual	1
Lanyard	1
Carrying case	1