

PH60-DE Replacement Electrode



Main Features

- Compatible with the Ph60-Z, PH60S-Z, PH60F-Z, and PC60-Z meters
- Highly sensitive bulb glass membrane with built-in temperature sensor
- Double junction design, suitable for dirty and demanding aqueous solutions, extending service life by up to 80%
- With a unique blue gel as inner buffer: no more air bubbles

Technical Parameters

pH	Measuring range	0 to 14
	Reference	Ag/AgCl
	Junction	Double ceramic
	Sensor type	Lithium glass
	Electrolyte	Gel KCl
Temperature	Measuring range	0 to 50 °C
	Temperature sensor	30 kΩ thermistor
Other	Body	Polycarbonate
	Dimensions	75×35×35 mm
	Weight	40 g