

PH60F-DE Replacement Electrode



Main Features

- Compatible with the PH60-Z, PH60S-Z, PH60F-Z and PC60-Z meters.
- Flat glass membrane with built-in temperature sensor.
- Double junction design, suitable for dirty and demanding aqueous solutions; extends service life by up to 80%.
- Ideal for surface pH measurement on skin, paper, fabric, and concrete, and for small-volume solutions.

Technical Parameters

pH	Measuring range	0 to 14
	Reference	Ag/AgCl
	Junction	Double ceramic
	Sensor type	Flat lithium glass
	Electrolyte	Gel KCl
Temperature	Measuring range	0 to 50 °C
	Temperature sensor	30 kΩ thermistor
Other	Body	Polycarbonate
	Dimensions	75×35×35 mm
	Weight	40 g