

PH60F-E Replacement Electrode



Main Features

- Compatible with the pH meters of the PH60 range
- pH determination on surfaces such as skin, paper, fabric, etc.
- Glass electrode with integrated temperature sensor
- Service life of up to 24 months with regular use

Technical Parameters

pH	Measuring range	0 to 14
	Reference	Ag/AgCl
	Junction	Ceramic
	Sensor type	Flat lithium glass
	Electrolyte	3M KCl
Temperature	Measuring range	0 to 50 °C
	Temperature sensor	30 kΩ thermistor
Other	Body	Polycarbonate
	Dimensions	75×35×35 mm
	Weight	40 g