

PH60S-DE Replacement Electrode



Main Features

- Compatible with the PH60-Z, PH60S-Z, PH60F-Z, and PC60-Z meters
- Pointed glass membrane with built-in temperature sensor
- Swiss sensor technology: long-life reference system with polymer electrolyte
- Ideal for pH measurement of cheese, meat, fruit, vegetables, bread, sushi rice, and other solid/semi-solid samples.

Technical Parameters

pH	Measuring range	0 to 14
	Sensor type	LabSen pointed glass membrane
	Electrolyte	Gel KCl
	Junction	Open
	Reference	Ag/AgCl
Temperature	Temperature sensor	30 k Ω thermistor
	Operating temp.	0~50°C (32 ~ 122°F)
Other	Body	Polycarbonate
	Dimensions	140×50×60 mm
	Weight	70 g