

PH700 Benchtop pH Meter



Main Features

- Includes the 201T-F 3-in-1 pH electrode for pH and temperature measurement, with automatic temperature compensation
- Automatic calibration mode with on-screen calibration guidance and self-diagnosis
- Modern digital technology delivers fast and stable readings
- Electrode holder and calibration solutions included
- Splash-proof (IP54 rating)



Technical Parameters

pH	Measuring range	0,00 to 14,00
	Resolution	0,1/0,01
	Accuracy	$\pm 0,01 \pm 1$ digit
	Calibration	Automatic 1- to 3-point calibration
	Electrode connector	BNC
	Temp.-Compensation	0 to 100 °C (automatic or manual)
mV	Measuring range	± 1999 mV
	Resolution	1 mV
	Accuracy	$\pm 0,1 \% \text{ F.S.} \pm 1$ digit
Temperature	Measuring range	0 to 100 °C
	Resolution	0,1 °C
	Accuracy	$\pm 0,5 \text{ °C} \pm 1$ digit
	Connector	RCA
Other	Memory capacity	50 groups
	Stored values	Number, reading, unit, temperature
	Power adapter	DC 9 V/ 300 mA
	IP rating	IP54
	Dimensions	24 × 23,5 × 10,3 cm (with case 30 × 25 × 8 cm)
	Weight	1 kg (with case 1,8 kg)

What's in the Kit

PH700 pH Meter	1
201T-F 3-in-1 pH electrode	1
pH4.00 Calibration solution (50ml)	1
pH7.00 Calibration solution (50ml)	1
pH10.01 Calibration solution (50ml)	1
Electrode holder	1
Power adapter	1
Manual	1