

SX610 pH Pocket Meter



Main Features

- Simple one-button operation
- Compact design with replaceable sensor; ideal for measurements in test tubes
- Smart functions: automatic shut-off, low-battery indicator
- Automatic stable-reading detection
- Complete kit incl. calibration solutions, storage solution, and batteries in a handy case



Technical Parameters

pH	Measuring range	0 to 14,0
	Resolution	0,1
	Accuracy	±0,1
	Calibration Points	Automatic 1-point calibration
	Temp.-Compensation	0 to 50 °C (automatic)
Temperature	Measuring range	0 to 50 °C
	Resolution	0,1 °C
	Accuracy	±0,5 °C
Other	Power supply	2 × CR2032 batteries
	Operating time	Up to 1000 Hours
	IP rating	IP67
	Dimensions	148×29×14 mm (with case 165×95×28 mm)
	Weight	49 g (with case 160 g)

What's in the Kit

SX610 pH Meter	1
pH7.00 Calibration solution	1
3M KCl Storage Solution	1
CR2032 battery	2
Manual	1
Case	1