

SX615 Replacement Electrode



Main Features

- pH determination of all kinds of liquids
- Glass electrode with integrated temperature sensor
- Slim electrode design (Ø10 mm), ideal for measurements in test tubes
- Service life of up to 24 months with regular use
- Compatible with the SX620 meter

Technical Parameters

pH	Measuring range	0,0 to 14,0
	Reference	Ag/AgCl
	Junction	Ceramic
	Sensor type	Lithium glass
	Electrolyte	3M KCl
Temperature	Measuring range	0 to 60°C
	Temperature sensor	30 kΩ thermistor
Other	Electrode	Ø10 mm
	Dimensions	100×30×15 mm
	Weight	16 g