

SX620 pH Pocket Meter



Main Features

- 0,01 pH accuracy
- Automatic 3-point calibration
- Automatic temperature compensation 0-60 °C
- A slim, replaceable probe suitable for test tubes
- Smart functions: switching between °F and °C, automatic shut-off
- Complete kit: 3 ready-to-use calibration solutions and 1 storage solution in a compact case



Technical Parameters

pH	Measuring range	-1 to 15,00
	Resolution	0,1
	Accuracy	±0,01
	Calibration Points	Automatic 1- to 3-point calibration
	Temp.-Compensation	0 to 50 °C (automatic)
Temperature	Measuring range	0 to 60 °C (automatic)
	Resolution	0,1 °C
	Accuracy	±0,5 °C
Other	Power supply	2 × CR2032 batteries
	Operating time	Up to 1000 Hours
	IP rating	IP67
	Dimensions	148×29×14 mm (with case 165×95×28 mm)
	Weight	49 g (with case 160 g)

What's in the Kit

SX620 pH Meter	1
pH4.00 Calibration solution	1
pH7.00 Calibration solution	1
pH10.00 Calibration solution	1
3M KCl Storage Solution	1
CR2032 battery	2
Manual	1
Case	1